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
- Naming the coronavirus disease (COVID-19) and the virus that causes it. (n.d.). Retrieved from [https://www.who.int/emergencies/diseases/novel-coronavirus-2019/technical-guidance/naming-the-coronavirus-disease-\(covid-2019\)-and-the-virus-that-causes-it](https://www.who.int/emergencies/diseases/novel-coronavirus-2019/technical-guidance/naming-the-coronavirus-disease-(covid-2019)-and-the-virus-that-causes-it).
- "Frequently Asked Questions - General Assembly of the United Nations." United Nations, www.un.org/pga/75/coronavirus/faqs/.
- World Health Organization. (2020). Advice on the use of masks in the community, during home care, and in health care settings in the context of COVID-19: interim guidance, 19 March 2020 (No. WHO/2019-nCoV/IPC_Masks/2020.2). World Health Organization.
- Healthcare Professionals: Frequently Asked Questions and Answers. (2020, March 22). Retrieved from <https://www.cdc.gov/coronavirus/2019ncov/hcp/faq.html>.

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PI-R0182C-ADC-FS Rev. D2.0
Date released: 2021-08-02
Vietnamese version

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Holds Certificate No:

FM 517416

and operates a Quality Management System which complies with the requirements of ISO 13485:2016 for the following scope:

Research, design and development, manufacturing and distribution of in vitro diagnostic reagents, controls, test kits and instruments for use in the areas of Autoimmune / Allergy, Cardiac Markers, Fertility, Infectious Diseases, Endocrinology, Oncology, and Veterinary.

For and on behalf of BSI:

Stewart Brain, Head of Compliance & Risk - Medical Devices

Original Registration Date: 2007-07-26

Latest Revision Date: 2019-05-10

Effective Date: 2019-05-16

Expiry Date: 2022-05-15

Page: 1 of 1

...making excellence a habit.™

This certificate remains the property of BSI and shall be returned immediately upon request.
An electronic certificate can be authenticated [online](http://www.bsigroup.com/ClientDirectory). Printed copies can be validated at www.bsigroup.com/ClientDirectory
To be read in conjunction with the scope above or the attached appendix.

Americas Headquarters: BSI Group America Inc., 12950 Worldgate Drive, Suite 800, Herndon, VA 20170-6007 USA
A Member of the BSI Group of Companies.

Certificate of CE-Registration

This is to certify that, in accordance with the *In Vitro* Diagnostic Medical Device Directive 98/79/EC, Medical Device Safety Service GmbH (MDSS) agrees to perform all duties and responsibilities as the Authorized Representative for:

**CTK Biotech, Inc.
13855 Stowe Dr.
Poway, CA 92064
USA**

as stipulated and demanded by the aforementioned Directive. The European Databank on Medical Devices (EUDAMED) is established as of May 1, 2011. The German Competent Authority is notified of the manufacturer's *in vitro* diagnostic medical devices and has allocated registration numbers shown in:

Annex A dated 2021-06-01

The Manufacturer has provided MDSS with the appropriate Declaration(s) of Conformity confirming that the *in vitro* diagnostic medical devices fulfill the applicable requirements of Directive 98/79/EC. In compliance with German law, a safety officer has been appointed for Germany.

2021-06-01

Dr. Philipp Hohenbrink
Director AR Services
MDSS GmbH

MDSS · Medical Device Safety Service · Schiffgraben 41 · 30175 Hannover, Germany

Annex A: 2021-06-01
Manufacturer: CTK Biotech, Inc.

REF	Device Names (notified)	Nomenclature (notified)	Code (notified)	Description (notified)	Class (notified)	EC Certificate No. & Expiry (notified)	German Registration Number
R0180C	OnSite COVID-19 IgG/IgM Rapid Test; Aria COVID-19 IgG/IgM Rapid Test; Answer COVID-19 IgG/IgM Rapid Test SURGITECH COVID-19 IgG/IgM Rapid Test	EDMA	15 70 90 08	Coronavirus - RT & POC	Other	N/A	DE/CA09/0170/C33/IVD/055-03
R0182C	OnSite COVID-19 Ag Rapid Test Aria COVID-19 Ag Rapid Test Answer COVID-19 Ag Rapid Test ad-bio COVID-19 Ag Rapid Test ImmuView COVID-19 Ag Rapid Test SURGITECH COVID-19 Ag Rapid Test	EDMA	15 70 90 08	Coronavirus- RT & POC	Other	N/A	DE/CA09/0170/C33/IVD/055-03
E0315	RecombiLISA CHIK IgM ELISA; Aria CHIK IgM ELISA	EDMA	15 04 80 25	Chikungunya Virus	Other	N/A	DE/CA09/0170/C33/IVD/019-03
E0320	RecombiLISA Zika IgM ELISA; Aria Zika IgM ELISA	EDMA	15 04 80 30	Zika Virus	Other	N/A	DE/CA09/0170/C33/IVD/031-03
R0188C	OnSite COVID-19 + Influenza A/B Ag Rapid Test Aria COVID-19 + Influenza A/B Ag Rapid Test	EDMA	15 04 80 90	Other Viral Antigen/Antibody Detection	Other	N/A	DE/CA09/0170/C33/IVD/072
E0170	RecombiLISA Chagas IgG ELISA; Aria Chagas IgG ELISA	EDMA	15 05 10 02	Chagas	Other	N/A	DE/CA09/0170/C33/IVD/018-04
E0121	RecombiLISA Leishmania Ab ELISA; Aria Leishmania Ab ELISA	EDMA	15 05 10 05	Leishmania	Other	N/A	DE/CA09/0170/C33/IVD/023-04
P0189	Aridia SARS-CoV-2 & Influenza A/B Real-Time PCR Test	EDMA	15 30 02 42	Influenza & Para Influenza / Respiratory Syncytial Virus (RSV) / SARS-CoV - Multiplex NA Reagents	Other	N/A	DE/CA09/0170/C33/IVD/079

U.S. FOOD & DRUG
ADMINISTRATION

CENTER FOR VETERINARY MEDICINE

CURRENT GOOD MANUFACTURING PRACTICE CERTIFICATE (CGMP)

CERTIFICATE NUMBER 6ZBH-GKQC	EXPORTING COUNTRY United States of America	IMPORTING COUNTRY No Country Applicable	EXPIRATION DATE February 09, 2023
FACILITY INFORMATION			
Manufacturer Name/Address CTK Biotech, Inc. 13855 Stowe Dr, Poway, California 92064 United States		FEI Number 3003952975	Date of Last Inspection 08/16/2018
PRODUCT INFORMATION (OPTIONAL)			
Rapid Test, ELISA, CLIA, RaFIA, and Instrumentation Devices			
ATTESTATION			
I am the Director, Division of Compliance, within the Center for Veterinary Medicine, Food and Drug Administration (FDA). In this capacity I issue Certificates to Foreign Government concerning the manufacture, preparation, and marketing of drugs in the United States and the Current Good Manufacturing Practice (cGMP) status of the plant that produces them. It is certified that the above manufacturing facility may manufacture, prepare, market and legally export products from the United States at this time. The plant is subject to the jurisdiction of the FDA and is subject to periodic inspections. The last such inspection showed that the plant, at this time, is in substantial compliance with Current Good Manufacturing Practice regulations as required by the Federal Food, Drug, and Cosmetic Act. This Certificate is not a substitute for a Certificate to Foreign Government that attests to the legality and exportability of specific animal drug products.			
NAME, TITLE Eric M. Nelson, Director, Division of Compliance, Center for Veterinary Medicine		DATE ISSUED February 09, 2021	
Signature: 			

This certificate expires 24 months from the date issued. To verify the authenticity of the information on this certificate, you may click www.access.fda.gov/fev/verifyCertificate to view details of the certificate as issued by FDA.

Page 1 of 1

EUROPEAN COMMISSION
DIRECTORATE-GENERAL FOR HEALTH AND FOOD SAFETY

Public health, country knowledge, crisis management
Health Security

EU health preparedness:

A common list of COVID-19 rapid antigen tests and a common standardised set of data to be included in COVID-19 test result certificates

Agreed by the Health Security Committee

This document was agreed by the HSC on 17 February 2021

Annex I

Common list of COVID-19 rapid antigen tests

A first update was agreed by the HSC on 10 May 2021; A second update was agreed by the HSC on 16 June 2021; A third update was agreed by the HSC on 7 July 2021; A fourth update was agreed by the HSC on 14 July 2021

IMPORTANT: A (interim) grace period of 8 weeks applies whenever updates are made to Annex I, the common list of COVID-19 rapid antigen tests

Annex II

Common standardised data set of to be included in COVID-19 test result certificates

An update to Annex II was agreed by the HSC on 19 March 2021

Manufacturer	RAT commercial name	CE marking	Clinical performance Data by manufacturer	Clinical performance Data used in MS	FIND evaluation studies	EU Member States using in practice	Other countries using in practice	Completed validation studies	MS currently validating	In JRC database (Device ID #) ¹¹
				DE: Positive evaluation by Paul-Ehrlich-Institut (sensitivity of 100% at <Ct25) + Manufacturer specificity: 100% NL: Independent field study, mainly symptomatic individuals, sensitivity Ct_≤30: 96.0%; specificity overall: 100%						
BIOSYNEX SA	BIOSYNEX COVID-19 Ag+ BSS	Yes	Clinical Sensitivity: 97.5 %	FR: Validation study data: 125 positive and 118 negative samples; sensitivity 96%, specificity: 99%		FR		FR		Yes (1494)
BIOTEKE CORPORATION (WUXI) CO., LTD	SARS-CoV-2 Antigen Test Kit (colloidal gold method)	Yes	96.49 % sensitivity 99.28 % specificity OP/NP swab	DE: Positive evaluation by Paul-Ehrlich-Institut (sensitivity of 95% at <Ct25) + Manufacturer specificity: 99.28%		DE ^[2]		DE ^[2]		Yes (2067)
BTNX Inc	Rapid Response COVID-19 Antigen Rapid Test	Yes	90.2% sensitivity 100% specificity NP swab, NP swab, OP swab	DE: 94.55% sensitivity, 100% specificity		AT, DE ^[2] , ES, SI		DE ^[2]		Yes (1236)
CerTest Biotec	CerTest SARS-CoV-2 Card test	Yes	92.9% sensitivity 99.6% specificity NP swab	ES: Ct ≤ 25, sensitivity: 94,0%; sensitivity for samples within the first 5 days after symptom onset: 84,8%		ES, PT, SI		DE ^[2] , ES		Yes (1173)
Core Technology Co., Ltd	Coretests COVID-19 Ag Test	Yes	98.1% sensitivity 99.6% specificity NP swab	DE: 98.1% sensitivity, 99.6% specificity		AT, DE ^[2] , RO		DE ^[2]		Yes (1919)
CTK Biotech, Inc	OnSite COVID-19 Ag Rapid Test	Yes	Clinical Sensitivity: 92.3 % Clinical Specificity: 100 % Nasal, NP swab	ES: 219 samples; Nasal swab - Clinical sensitivity 86% (90%: Ct <30) Specificity: 100% (Method B) DK: 107 samples; Nasal swab - clinical sensitivity 86%; (from asymptomatic and mild symptomatic individuals), Clinical specificity: 100%	To start	DK		DK, ES		Yes (1581)
DDS DIAGNOSTIC	Test Rapid Covid-19 Antigen (tampon nazofaringian)	Yes	98.77% sensitivity 99.03% specificity Nasal swab	RO: Meets the minimum performance requirements.		RO		RO China	RO	Yes (1225)

25.06.2021

Comparative evaluation of the sensitivities of SARS-CoV-2 antigen rapid tests

Aim

Comparison of different antigen rapid tests with using identical sample material

Material

Pools from nasopharyngeal and oropharyngeal swabs.

Dry swabs were included in PBS; moist swabs were already included in the transport media of various compositions. Pools are random mixtures obtained from up to 10 samples of comparable CT values diluted 1:10 in negative samples in PBS. The CT values of a pool were determined by means of different PCR assays, and the putative number of RNA copies calculated with the aid of the INSTAND standards. In the case of the PCRs used, a CT value of 25 corresponds to around 10^6 RNA copies/mL. 18 samples each were analysed with $CT < 25$, 23 samples with CT between 25 and 30, and 9 samples with $CT > 30$. The replication of the virus in cell culture was determined as a possible correlate for infectiousness as another characteristic of the samples.

Method

The pools were aliquoted, frozen, shipped, and thawed for evaluation of the tests. For each test, 50 µL of the pool were analysed using the components of the test provided, e.g. swabs. Laboratories participating in the comparative evaluation included the Robert Koch-Institut, the Paul-Ehrlich-Institut, the reference laboratory for coronaviruses (Charité), and the Institute for Microbiology of the German Army (Bundeswehr).

Summary

This comparative evaluation of a large number of SARS-CoV-2 rapid antigen tests (point of care tests; POCT) of different designs and manufacturers with the same sample set allows an overview of the current state of art regarding sensitivity. The results do not allow any conclusions regarding specificity of the tests.

Those POCTs which have up to now been included in the evaluation and have been assessed as reflecting the current state of the art are listed in the table below. Other tests, which were assessed as not reflecting the state of the art were deleted from the list of the Federal Institute for Drugs and Medical Devices (Bundesinstitut für Arzneimittel und Medizinprodukte, BfArM). This comparative evaluation is constantly continued, and the table is amended accordingly.

You should be aware that this comparative evaluation can only cover a random sample of the SARS-CoV-2 rapid antigen tests listed by the BfArM, thus eligible for refunding, and that few other products could not (yet) be taken into account, despite the interests on the part of the manufacturers/distributors.

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**Certificate of Marketability
of Medical Devices**

Confirmation of marketability of medical devices according to § 34 (1) of the Medical Devices Law (MPG) of 07. August 2002 (BGBl. I S. 3146) in the version valid until 25. May 2021.

For presentation to the competent authorities or bodies of **Viet Nam**

It is certified that the following medical devices can be marketed without restriction within

- Germany
- the member states of the European Union and
- the other states having a contractual agreement with the European Economic Area.

Products:

see Annex A

Manufacturer:

CTK Biotech, Inc.
13855 Stowe Dr.
Poway, CA 92064
United States of America

Authorized Representative:

Medical Device Safety Service GmbH (MDSS)
Schiffgraben 41
30175 Hannover
Germany

It is also certified that with the affixing of the CE-mark, the manufacturer follows the prescribed conformity assessment procedure and ensures and declares that the medical devices meet the essential requirements of the Council Directives 98/79/EC of 27. October 1998 in the present valid version.

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Anlage A / Annex A

Product	Registration	EDMA
OnSite COVID-19 IgG/IgM Rapid Test: R0180C	DE/CA09/0170/C33/IVD/055-04	15 70 90 08
OnSite COVID-19 Ag Rapid Test: R0182C		
OnSite COVID-19 + Influenza A/B Ag Rapid Test: R0188C	DE/CA09/0170/C33/IVD/072	15 04 80 90

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Im Auftrage

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Australian Government

Department of Health
Therapeutic Goods Administration

Australian Register of Therapeutic Goods Certificate

Issued to

MD Solutions Australasia Pty Ltd

for approval to supply

Severe acute respiratory syndrome-associated coronavirus IVDs

ARTG Identifier	332961
ARTG Start Date	31/03/2020
Product Category	Medical Device Included - IVD Class 3
GMDN	CT772
GMDN Term	Severe acute respiratory syndrome-associated coronavirus IVDs
Intended Purpose	The OnSite/Aria COVID-19 IgG/IgM Rapid Test is a lateral flow immunoassay for the detection of anti-SARS-CoV-2 IgG and IgM antibodies in human serum, plasma or whole blood. It is intended to be used by healthcare professionals as an aid in the diagnosis of infection with SARS-CoV-2 coronavirus, which causes COVID-19 disease. Any interpretation or use of this preliminary test result must also rely on other clinical findings as well as on the professional judgment of healthcare providers. Alternative test method(s) should be considered to confirm the test result obtained by this device.

Manufacturer Details	Address	Certificate number(s)
CTK Biotech Inc	13855 Stowe Drive Poway , California , 9 2064 United States Of America	DV-2020-MC-03953-1

ARTG Standard Conditions

The above Medical Device Included - IVD Class 3 has been entered on the Register subject to the following conditions:

- The inclusion of the kind of device in the ARTG is subject to compliance with all conditions placed or imposed on the ARTG entry. Refer Part 4-5, Division 2 (Conditions) of the Therapeutic Goods Act 1989 and Part 5, Division 5.2 (Conditions) of the Therapeutic Goods (Medical Devices) Regulations 2002 for relevant information.
- Breaching conditions of the inclusion related to the device of the kind may lead to suspension or cancellation of the ARTG entry; may be a criminal offence; and civil penalties may apply.

Products Covered by This Entry

1. Severe acute respiratory syndrome-associated coronavirus IVDs

This entry: contains System(s)/Procedure Pack(s)

IVD Information

Name	Category Description
OnSite COVID-19 Ag Rapid Test	Point of care testing
Aria COVID-19 IgG/IgM Rapid Test	Point of care testing
Aria COVID-19	Point of care testing

Name	Category Description
Ag Rapid Test	
OnSite COVID-19 IgG/IgM Rapid Test	Point of care testing

Product Specific Conditions

No specific conditions have been recorded against this entry.

Therapeutic Goods Administration
PO Box 100, Woden ACT 2606 Australia
Phone: 1800 020 653
Email: info@tga.gov.au

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Att. Shauna Parker
CTK Biotech A/S
13855 Stowe Dr, Poway
CA 92064
USA

26 March 2021
Case no. 2021033589
T +45 44 88 95 95
E med-udstyr@dkma.dk

The Danish Medicines Agency hereby authorises CTK Biotech A/S to place on the market/put into service the OnSite COVID-19 Ag Rapid Test in Denmark.

The authorisation is conditional with the following terms and conditions.

Terms and conditions:

- This authorisation is valid until 01.10.2021, or until the date when the device is **CE marked**, if this occurs before 01.10.2021. The manufacturer must inform the Danish Medicines Agency when the CE-mark has been acquired.
- The authorisation is limited to marketing/putting into service under supervision in Danish schools and Danish educational institutions.
- If the CE marking has not been achieved by 01.07.2021, the manufacturer shall provide a status for the CE marking process.
- **CTK Biotech A/S** (manufacturer) shall establish a vigilance system, that ensures the collection and reporting to **the Danish Medicines Agency** of any incident, which occurs during the authorisation period.
- The manufacturer must ensure traceability of products through, for example, batch / LOT numbers.
- The information on the labelling and in the instruction for use (IFU) for the "OnSite COVID-19 Ag Rapid Test" shall be in Danish, before distribution to the end-user.
- The labelling must clearly state that the products are placed on the market under this authorisation as a derogation from the conformity assessment procedures.
- The test must be accompanied by a detailed IFU, intended specifically for lay people and the intended user group (Danish educational institutions under supervision).
- Technical documentation for the product must be kept for a minimum of 5 years.
- A status shall be provided 01.07.2021 and by the end of this authorisation period of how many devices have been placed on the market under this authorisation.

The authorisation covers the following product:

Product: OnSite COVID-19 Ag Rapid Test
Catalogue/ Model number(s): R0182C
In Vitro Medical Device: General IVD, non-CE-marked for self-testing
Scope: lateral flow immunoassay for the qualitative detection of SARS-CoV-2 nucleocapsid (NP) antigens in nasal swab specimens

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Manufacturer of the medical device: CTK Biotech A/S
13855 Stowe Dr, Poway
CA 92064
USA

Authorised representative for manufacturer: MDSS -Medical Device Safety Service
Schiffgraben 41
30175 Hannover
Germany

Summary of the application of authorisation:

On the 23th of March 2021 CTK Biotech A/S has applied for an authorisation to place on the market/put into service a non-CE-marked in vitro diagnostic medical device for self-testing by way of derogation from the conformity assessment procedures in the Danish executive order for in vitro medical devices.

The manufacturer has in their application informed the Danish Medicines Agency that OnSite COVID-19 Ag Rapid Test is currently undergoing CE marking for self-testing purposes by the Notified Body TÜV Rheinland LGA Products GmbH, No. 0197, Nürnberg, Germany.

The application states that the tests are intended for use in schools and educational institutions for self-testing under supervision. The manufacturer will supply the products with a Danish IFU.

Reason:

The Danish Medicines Agency can according to section 6, subsection 11 in the Executive Order 1269 of 12. December 2005 on in vitro diagnostic medical devices, under certain circumstances grant an authorisation to market and put into service an in vitro diagnostic medical device, for which the procedures referred to in subsections 1 to 4 have not been carried out, for the protection of health.

The Danish Medicines Agency only considers applications if there are clear health reasons and if there are no available alternatives on the market.

In this case the Danish Health Authority has stated, that in the present circumstances related to the Covid-19 situation in Denmark, there is a need for students in Danish schools and educational institutions to be tested on a regular basis by using antigen test for self-testing. This is in order to maintain epidemic control and to contribute to improvement of the mental health of the students by letting them return to schools.

The Danish Health Authority has further stated, that it is not possible to organise testing of all students on a regular basis during a longer period, solely by using the present available national test system, where only tests that are CE marked for professional use are used.

After a review of OnSite COVID-19 Ag Rapid Test, Statens Serum Institut (SSI), which is the Danish agency for disease preparedness and surveillance, has stated that based on the data delivered by the manufacturer, the product is sufficient for the purpose of maintaining epidemic control.

It is therefore the Danish Medicines Agency's opinion, that there is a necessity for the protection of health of Danish students and that there are no available alternatives on the market at the present time.

Based on the statement from SSI, it is furthermore the Danish Medicines Agency's opinion, that the product OnSite COVID-19 Ag Rapid Test is sufficient for the purpose of maintaining epidemic control schools and educational institutions for self-testing under supervision.

The Danish Medicines Agency therefore authorises CTK Biotech A/S to place on the market/put into service the OnSite COVID-19 Ag Rapid Test in Denmark for the protection of the public health.

Legal basis

Medical devices are regulated in the EU by the Council Directive 98/79/EEC on in vitro diagnostic medical devices. The Directive is implemented in the Danish Executive Order No. 1269 of 12 December 2005 on In Vitro Diagnostic Medical Devices.

This decision is made in accordance with section 6, subsection 11 of the Executive Order on in vitro diagnostic medical devices, that states, that The Danish Medicines Agency, by way of derogation from subsections 1 to 4, may authorise, on duly justified request, the placing on the market and putting into service, of individual devices for which the procedures referred to in paragraphs 1 to 4 have not been carried out and the use of which is in the interest of protection of health.

Appeals procedure

This decision may be appealed to the Danish Ministry of Health, Holbergsgade 6, 1057 Copenhagen K, sum@sum.dk.

Sincerely

Thomas Wejs Møller
Head of Unit

Jeppe Larsen
Head of Section

MINISTRY OF HEALTH
PHARMACY AND POISONS BOARD

Telegram: "MINHEALTH" Nairobi
Telephone: 020-2716905/6, 020-3562107
Cellphone: 0733-884411/0720 608811
Fax: 2713409
Email: admin@pharmacyboardkenya.org
Website: www.pharmacyboardkenya.org

Pharmacy & Poisons Board Hse
Along Lenana Road
P. O. Box 27663-00506
NAIROBI

When replying please quote our ref No:
PPB/PER/VOL.1/LET/031/21

6th April 2021

MAH (Manufacturer/Supplier)

CTK Biotech Inc.
13855 Stowe Drive, Poway
CA 92064 USA

Through'

Local Technical Representative

Pharmadex International Ltd
P.O. Box 53476 - 00200
Nairobi, Kenya

Dear Sir/Madam

**RE: APPLICATION FOR IN-VITRO DIAGNOSTIC MEDICAL DEVICE
EMERGENCY USE AUTHORIZATION OF CTK Biotech Onsite COVID-19
Rapid Antigen Test**

Reference is made to the above subject.

The Pharmacy and Poisons Act (Cap 244) as amended by the Health Laws (Amendment) Act, 2019, provides for the regulatory oversight of Medical Devices including In-Vitro Diagnostics (IVDs) in Kenya by the Pharmacy and Poisons Board "the Board". Please note that the Expert committee on medicals devices & IVDs of the Board is tasked with reviewing applications for serological test kits (Covid-19). The Board's policy requires that all IVDs meet certain conditions prior to approval; including submission of product dossier and validation protocols and reports that evidence threshold of key performance parameters. The validating laboratory must be an accredited laboratory, such as, Kenya Medical Research Institute that is recognized by the Board.

The Committee upon review of your application in a meeting held on 31st March, 2021 passed a verdict of Recommended for Emergency Use Authorization. *Note that use of this device for diagnosis of Covid-19 or any other unauthorized purpose is strictly prohibited and may lead to revocation of this license and or criminal liability to the Applicant (**Local Technical Representative**) and any other concerned party.*

The EUA will be revoked due to any reason that may necessitate the regulatory action or once the Covid-19 Pandemic is declared over by the ministry of Health, or the Head of State whichever one comes first; In the later scenario, you will be required to submit the clinical information for full registration of your product in the Kenyan Market.

Please see overleaf the details of the authorized product.


Yours faithfully

for

Dr. Ronald M. Inyangala

Deputy Director, Product Evaluation and Registration
For: CHIEF EXECUTIVE OFFICER

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ANNAR DIAGNOSTICA IMPORT S.A.S. con domicilio en BOGOTÁ - D.C.